DVIA-T Series

Tabletop Active Vibration Isolation Platform

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Contents


1. Components	4
2. Installation Procedure	4 - 5
3. How to Active the System	6 - 7
4. How to Turn Off the System	7
5. How to Lock the System	8
7. Troubleshooting	9
8. Warranty	10

1. Components


- 1) Tabletop Active Vibration Isolation Platform
- 2) AC Power Cable
- 3) RS232 Cable

2. Installation Procedure


- 1) Place DVIA-T on a rigid table (A table must be rigid as possible and should not be wobbly).
- 2) Adjust leveling foot until DVIA-T is firmly placed on a table.

- 3) Place an instrument on DVIA-T. The center of the instrument mass should be located in the center of DVIA-T.

- 4) Make sure that a power switch is off, then connect an AC Power Cable to Power Inlet at the back side of DVIA-T and plug the AC Power Cable to an electrical outlet.

On

Off

3. How to Active the System



- 1) Press a Power Switch.
- 2) Press an Auto Leveling Button and release.


During the auto leveling mode, the top plate goes up and down automatically to adjust tis level.

Please avoid any physical impacts on the top plate during the auto leveling.

- 3) Check level meters on LCD display. The following figures indicate positions of Z1, Z2, Z3 and Z4. If one of them relatively higher than others, the center of mass the instrument is not evenly distributed.

Normal Level Meter

Abnormal Level Meter

When the level meters are abnormally shown, relocate the instrument in the center on DVIA-T.

- 4) When the auto leveling is completed, the auto leveling button's light is off and the LCD display is shown as the figure below.

- 5) Press an active isolation button and release.

4. How to Turn Off

- 1) Before Turning off DVIA-T, disable the current active vibration isolation mode first.
- 2) Press an active vibration isolation button to disable the active vibration isolation mode
- 3) Press a power switch button to turn off the system.

5. How to Lock the System


- 1) Unload an instrument from DVIA-T
- 2) Press an active isolation button and auto leveling button at the same time and then release.

- 3) DVIA-T is internally locked for transportation.
- 4) To activate the system, follow steps in **3. How to Activate the Systems.**


6. Troubleshooting

- 1) Problem: After tuning on DVIA-T, motors do operate with the following message shown on the LCD display when an auto leveling button is pressed.


Solution: This problem happens when the processor cannot read the input of the position sensor properly. This problem can be solved by restarting DVIA-T.

- 2) Problem: When pressing an active isolation button, vibration signals are shown abnormally on the LCD display like the following figures.

Normal Vibration Level

Abnormal Vibration Level

Solution: This problem happens when the processor cannot read the input of the position sensor properly. This problem can be solved by restarting DVIA-T.

Please contact DAEIL SYSTEMS, if you experience any issues with our products.

Email: internationalsales@daeilsys.com

DAEIL SYSTEMS LIMITED WARRANTY

DAEIL SYSTEMS warrants any defects in material and workmanship under normal use for one year from the date of original retail purchase by the original end-user purchaser. During this warranty period, DAEIL SYSTEMS will only repair or replace DAEIL SYSTEMS products at no cost only if products are not subject to abuse, misuse, neglect, carelessness, accident, alteration, fire, repair or other external factors. To employ DAEIL SYSTEMS warranty service, contact DAEIL SYSTEMS to get the best course of action. Replaced or repaired DAEIL SYSTEMS products are warranted for the remainder of the original warranty period.

Limitation of Warranty

TO THE EXTENT PERMITTED BY LAW, THE WARRANTY AND REMEDIES SET FORTH ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES AND REMEDIES, AND WE SPECIFICALLY DISCLAIM ALL STATUTORY OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND AGAINST HIDDEN OR LATENT DEFECTS. IF WE CANNOT LAWFULLY DISCLAIM STATUTORY OR IMPLIED WARRANTIES, THEN TO THE EXTENT PERMITTED BY LAW, ALL SUCH WARRANTIES SHALL BE LIMITED IN DURATION TO THE DURATION OF THIS EXPRESS LIMITED WARRANTY AND TO REPAIR OR REPLACEMENT SERVICE AS DETERMINED BY DAEIL SYSTEMS IN ITS SOLE DISCRETION.

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Manufacturer : DAEIL SYSTEMS CO.,LTD