


NOTES:


1. CONSTRUCTION:
Stainless Steel (400 Series) Top and Bottom Plates, Plated Steel Honeycomb Core, High Pressure Laminate Core Sidewalls
2. CORE CELL AREA: $\leq 3 \text{ cm}^2$
3. WEIGHT: 34.00 kg
4. ROHS: COMPLIANT



SPECIFICATIONS SUBJECT TO CHANGE WITHOUT NOTICE
DIMENSIONS ARE FOR REFERENCE ONLY

**FOR INFORMATION ONLY:
DO NOT MANUFACTURE
PARTS TO THIS DRAWING**

EO® **Edmund Optics**®

THIRD ANGLE PROJECTION 		TITLE	600 x 1200mm, Metric Lab Table	
ALL DIMS IN	mm	DWG NO	55252	SHEET 1 OF 1