

High Performance Vibration Isolation System



Stock #87-226

⊖ 1 ⊕ **\$3,854.00**

ADD TO CART

Qty 1-4

\$3,854.00

Qty 5+

\$3,662.00

Volume Pricing

[Request Quote](#)

Type:
Universal

Thickness (inches):
3.0

Length (inches):
20.0

Width (inches):
18.0

Construction:

Electropolished stainless steel sides and base
Textured powder coat carbon steel top plate

Feed Pressure (psi):

80 (nitrogen or air)

Height (inches):

3.0

Load Capacity (kg):

67 (150 lb)

Weight (kg):

22 (50 lb)

Regulatory Compliance

RoHS 2015:

[Compliant](#)

Certificate of Conformance:

[View](#)
