




25.0mm  
[0.98in]  
TYP



5.0mm  
[0.20in]

5.0mm  
[0.20in]



4 x M6 x 1.0 THRU SKIN

FOR INFORMATION ONLY  
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION		
NAME	DATE	
DRAWN GRM	15/MAY/13	
APPROVAL CJL	15/MAY/13	
COPYRIGHT © 2013 BY THORLABS		
VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUND OFF ERRORS		

<b>THORLABS</b> www.thorlabs.com	
NEXUS BREADBOARD	
MATERIAL	N/A
ITEM #	B90150B
APPROX WEIGHT	155 kg
REV	B